



PRODUCT PERFORMANCE (STRUCTURAL COMPOSITION):

BFE-95% - 25gsm

This mask is composed of non-woven 3-PLY – melt-blown - spun polypropylene whose soft layers are specifically engineered to comply with critical particulate filtration standard commonly recommended by health authorities for occupational or general use during times pandemic or environmental crisis.

SCOPE OF APPLICATION:

This mask is suitable for all types of personal use, while covering the user's mouth, nose and jaw. This mask will function to provide a specific physical barrier to prevent dust, as well as certain pathogenics, microorganisms, particulate matter, etc. from directly penetrating the mask itself.

INSTRUCTION FOR USE:

Unfold the mask with the metal nose clip facing upwards, then attach the ear loops around the ear and adjust to fit into the nose clip.

PRECAUTIONS:

Check the packaging before use. Damaged packaging VOIDS mask for use. This mask is for ONE-TIME use only. If you have any discomfort such as breathing difficulty- immediately remove and discontinue use of mask.

DISPOSAL:

This mask is biodegradable, however its use should dictate disposal options. Masks exposed or thought to be exposed to viral or bacterial contaminants should be considered and treated accordingly.

STORAGE CONDITIONS:

This product should be stored in a non-corrosive gas, cool, dry, well-ventilated environment, and away from flammable materials.

DIMENSIONS:

6.89in x 3.74in (17.5cm x 9.5cm)

PRODUCTION DATE, PRODUCTION BATCH NUMBER, VALIDITY DATE:

See packaging

CERTIFICATION:

FDA and CE

**DISPOSABLE
FOR ONE USE ONLY**



USAGE: As shown below

